

## Project NGS2468 data delivery, 09 Jun 2021:

Samples: 61  
Sequencing amount: 91.5 million read pairs  
Sequencing type: 150bp paired-end read (Illumina NextSeq 500/550 v2)  
Restriction enzyme(s): MslI

### Delivery contents:

- 'RAW': raw sequencing data after basecalling and demultiplexing in compressed FASTQ format
- 'AdapterClipped': compressed FASTQ files containing sequencing adapter clipped reads
- 'RE\_processed': compressed FASTQ files containing restriction enzyme site filtered reads

FastQC v0.11.9 [1] reports, containing read quality metrics, are stored along with the FASTQ files.

### Data analysis:

- Demultiplexing of all libraries for each sequencing lane using the Illumina bcl2fastq v2.20 software [2] (folder 'RAW', 'Group' subfolders):
  - 1 or 2 mismatches or Ns were allowed in the barcode read when the barcode distances between all libraries on the lane allowed for it
- Demultiplexing of library groups into samples according to their inline barcodes and verification of restriction site (folder 'RAW', 'Sample' subfolders):
  - no mismatches or Ns were allowed in the inline barcodes, but Ns were allowed in the restriction site
- Clipping of sequencing adapter remnants from all reads (folder 'AdapterClipped'):
  - reads with final length < 20 bases were discarded
- Restriction enzyme site filtering of read 5' ends (folder 'RE\_processed'):
  - reads with 5' ends not matching the restriction enzyme site are discarded
- Creation of FastQC reports for all FASTQ files
- Generation of read\_counts.xlsx, containing all read counts for all samples at a glance

If you have any questions related to your data or some steps of the data analysis, do not hesitate to contact us directly:

Email: [bioinformatics@lgcgroup.com](mailto:bioinformatics@lgcgroup.com)

Tel.: +49 (0)30 5304 2202

## References

- [1] Simon Andrews. *FastQC A Quality Control tool for High Throughput Sequence Data*. URL: <http://www.bioinformatics.babraham.ac.uk/projects/fastqc/>.
- [2] Illumina. *bcl2fastq2 Conversion Software*. URL: [https://support.illumina.com/sequencing/sequencing\\_software/bcl2fastq-conversion-software.html](https://support.illumina.com/sequencing/sequencing_software/bcl2fastq-conversion-software.html).